



**FRESENIUS
KABI**

caring for life

Freka[®] ENFit[®] Enteral Feeding Tube

PRODUCT RANGE





Contents



GASTRIC FEEDING TUBES	02-07	ACCESSORY ITEMS	10-13
> Freka CH/FR 8 Nasogastric Feeding Tube 80cm	02	> Freka Endo Introducer (280cm)	07
> Freka CH/FR 8 Nasogastric Feeding Tube 120cm	02	> Freka Click adaptor for Freka PEG Intestinal	10
> Freka CH/FR 10 Nasogastric Feeding Tube 120cm	03	> Freka Repair Set for Freka Surgical Jejunostomy, ENFit	10
> Freka CH/FR 12 Nasogastric Feeding Tube 120cm	03	> Freka Y Connector for Freka Intestinal	10
> Freka GastroTube CH/FR 15	04	> Freka Stoma Measuring Device	10
> Freka Button CH/FR 15 Gastrostomy	04	> Freka PEG Connector CH/FR 9 ENFit	11
> Freka Belly Button Gastrostomy	05	> Freka PEG Connector CH/FR 15 ENFit	11
> Freka PEG Standard Gastric CH/FR 9	05	> Freka PEG Connector CH/FR 20 ENFit	11
> Freka PEG Standard Gastric CH/FR 15	06	> Freka Cap ENFit	11
> Freka PEG Standard Gastric CH/FR 20	06	> Fixation plate for PEG, CH/FR 9	12
> Freka PEXACT CH/FR 15	07	> Fixation plate for PEG, CH/FR 15	12
JEJUNAL FEEDING TUBES	07-09	> Fixation plate for PEG, CH/FR 20	12
> Freka Endolumina Nasojejunal Feeding Tube	07	> Repair set for Freka PEG, CH/FR 9	12
> Freka Easy In Nasojejunal Feeding Tube	08	> Repair set for Freka PEG, CH/FR 15	12
> Freka Trelumina Nasojejunal Feeding Tube	08	> Repair set for Freka PEG, CH/FR 20	12
> Freka Surgical Jejunostomy Tube	09	> Freka Extension set (30cm)	13
> Freka PEG Universal Intestinal Feeding Tube	09	> Freka Extension set (100cm)	13
		> Cone adaptor	13
		> EasyBag Bolus adaptor	13
		ORDERING INFORMATION	14



All Freka feeding tubes and accessory items are ENFit compliant.



(ENFit products in this brochure are compliant with ENFit ISO 80369-3)


GASTRIC FEEDING TUBES

	<p>Description</p> <p>Product code: 7981839</p> <ul style="list-style-type: none"> > CH/FR 8, 80cm > ENFit connector > External diameter: 2.8mm > Internal diameter: 2.1mm 	
<p>NASOGASTRIC OR NASOJEJUNAL FEEDING TUBE</p> <p>FREKA</p> <p>NASOGASTRIC FEEDING TUBE</p>	<p>Key Features</p> <ul style="list-style-type: none"> > Tissue compatible tube made from polyurethane with ENFit Connector and integral cap > Spiral guide wire with connection for the injection of radiological contrast medium or lubricating fluid > Flexible polyurethane olive tip with 2 openings at the distal tube end > 3 radio-opaque strips offering enhanced visibility > 1cm interval markings along length of tube > Three-piece nasal adhesive tape set for secure tube fixation > Tube marked to confirm size, manufacturer and article code > Outlet holes in tip to help prevent kinking > Latex, PVC, DEHP free 	

	<p>Description</p> <p>Product code: 7981842</p> <ul style="list-style-type: none"> > CH/FR 10, 120cm > ENFit connector > External diameter: 3.4mm > Internal diameter: 2.2mm 	
<p>NASOGASTRIC OR NASOJEJUNAL FEEDING TUBE</p> <p>FREKA</p> <p>NASOGASTRIC FEEDING TUBE</p>	<p>Key Features</p> <ul style="list-style-type: none"> > Tissue compatible tube made from polyurethane > Spiral guide wire with connection for the injection of radiological contrast medium or lubricating fluid > Flexible polyurethane olive tip with 2 openings at the distal tube end > 3 radio-opaque strips offering enhanced visibility > 1cm interval markings along tube length > Three-piece nasal fixation tape set for secure tube fixation > Tube marked to confirm size, manufacturer and article code > Outlet holes in tip to help prevent kinking > Latex, PVC, DEHP free 	


	<p>Description</p> <p>Product code: 7981840</p> <ul style="list-style-type: none"> > CH/FR 8, 120cm > ENFit connector > External diameter: 2.8mm > Internal diameter: 2.1mm 	
<p>NASOGASTRIC OR NASOJEJUNAL FEEDING TUBE</p> <p>FREKA</p> <p>NASOGASTRIC FEEDING TUBE</p>	<p>Key Features</p> <ul style="list-style-type: none"> > Tissue compatible tube made from polyurethane with ENFit Connector and integral cap > Spiral guide wire with connection for the injection of radiological contrast medium or lubricating fluid > Flexible polyurethane olive tip with 2 openings at the distal tube end > 3 radio-opaque strips offering enhanced visibility > 1cm interval markings along length of tube > Three-piece nasal adhesive tape set for secure tube fixation > Tube marked to confirm size, manufacturer and article code > Outlet holes in tip to help prevent kinking > Latex, PVC, DEHP free 	

	<p>Description</p> <p>Product code: 7981843</p> <ul style="list-style-type: none"> > CH/FR 12, 120cm > ENFit connector > External diameter: 3.9mm > Internal diameter: 2.9mm 	
<p>NASOGASTRIC OR NASOJEJUNAL FEEDING TUBE</p> <p>FREKA</p> <p>NASOGASTRIC FEEDING TUBE</p>	<p>Key Features</p> <ul style="list-style-type: none"> > Tissue compatible tube made from polyurethane > Spiral guide wire with connection for the injection of radiological contrast medium or lubricating fluid > Flexible polyurethane olive tip with 2 openings at the distal tube end > 3 radio-opaque strips offering enhanced visibility > 1cm interval markings along tube length > Three-piece nasal fixation tape set for secure tube fixation > Tube marked to confirm size, manufacturer and article code > Outlet holes in tip to help prevent kinking > Latex, PVC, DEHP free 	



Description Product code: 7755648

- > Gastric CH/FR 15
- > ENFit Connector
- > 13cm shaft length
- > External diameter: 5mm
- > Internal diameter: 3.6mm



REPLACEMENT GASTROSTOMY


FREKA

GASTROTUBE BALLOON RETAINED REPLACEMENT GASTROSTOMY TUBE

For long-term intragastric feeding


Key Features

- > Latex free silicone tube - balloon retention
- > ENFit Connector
- > Rounded distal opening
- > Long life span balloon
- > Introducer device
- > 10ml luer syringe (for tube inflation only)
- > Tube clamp
- > Radio-opaque ring and strip
- > Step connector
- > Latex, PVC, DEHP free



Description Freka Belly Button Gastrostomy

- > CH/FR 12, 14, 16, 18, 20, 24
- > External diameter: 12(4mm), 14(4.7mm), 16(5.3mm), 18(6mm), 20(6.7mm), 24(8mm)
- > Internal diameter: 12(2mm), 14(2.4mm), 16(3mm), 18(3.6mm), 20(4mm), 24(5mm)



REPLACEMENT GASTROSTOMY LOW PROFILE

FREKA

BELLY BUTTON BALLOON GASTROSTOMY RETAINED


For medium to long-term intragastric feeding using a low-profile feeding device

Freka Belly Button	CH/FR size/shaft length	CH/FR 12	CH/FR 14	CH/FR 16	CH/FR 18*	CH/FR 20*	CH/FR 24*
Containing:	1.0 cm	7751612	7751622	7751632	7751642	7751652	7751662
1 Freka Belly Button	1.5 cm	7751613	7751623	7751633	7751643	7751653	7751663
1 Introducer (simple)	2.0 cm	7751614	7751624	7751634	7751644	7751654	7751664
1 Syringe for balloon filling	2.5 cm	7751615	7751625	7751635	7751645	7751655	7751665
10 ml Ampoule Aqua dest.	3.0 cm	7751616	7751626	7751636	7751646	7751656	7751666
1 Lubricant	3.5 cm	7751617	7751627	7751637	7751647	7751657	7751667
2 Gauze Pads	4.0 cm	7751618	7751628	7751638	7751648	7751658	7751668
	4.5 cm	7751619	7751629	7751639	7751649	7751659	7751669
	5.0 cm	7751620	7751630	7751640	7751650	7751660	7751670

Key Features

- > Medical grade silicone
- > Flat distal end
- > Valve for balloon filling
- > Flat and small button head
- > Anti-reflux valve
- > Button introducer device
- > Syringe plus 10ml ampoule for balloon filling
- > Latex, PVC, DEHP free

* Freka Belly Button sizes CH/FR 18 CH/FR 20 CH/FR 24 may have longer delivery times.




Description

Freka Button Gastrostomy Single Set

- > CH/FR 15
- > External diameter: 4.5mm
- > Internal diameter: 3.6mm
- > Shaft lengths (Product codes):

1.1cm (7751241) • 1.4cm (7751251) • 1.7cm (7751261)
2.0cm (7751271) • 2.5cm (7751281) • 3.0cm (7751291)
3.5cm (7751301) • 4.0cm (7751311) • 4.5cm (7751321)



REPLACEMENT GASTROSTOMY LOW PROFILE

FREKA

BUTTON BALLOON RETAINED

For long-term intragastric feeding using a low profile feeding device


Key Features

- > Medical grade silicone button with single bonded seam for durability¹
- > Long life span balloon
- > Horizontal and vertical radio-opaque strips
- > Flat open distal end
- > Anti-reflux valve
- > Valve for balloon filling
- > Retention balloon with 7.5ml fluid filled capacity
- > Latex, PVC, DEHP free

SINGLE SET


Single button set contains:

- > CH/FR 15 button with specific shaft length
- > 10ml luer syringe (for tube inflation only)
- > Button introducer device



Description Product code: 7755642

- > CH/FR 9
- > ENFit connector
- > Fine lumen, gastric
- > External diameter: 2.9mm
- > Length: 30cm
- > Internal diameter: 1.9mm



PERCUTANEOUS ENDOSCOPIC GASTROSTOMY USING THE PULL TECHNIQUE




FREKA




PEG STANDARD GASTRIC CH/FR 9

For long-term intragastric feeding




Key Features




- > Polyurethane tube with radio-opaque strip and silicone disc
- > Radio-opaque retention disc made from silicone
- > Puncture cannula with safety (air) valve
- > Introducer cone with double thread
- > Scalpel
- > External colour and size coded fixation plate (silicone) for securing tube (white/yellow)
- > ENFit connector
- > Interval markings every 2cm (from 2cm - 16cm)
- > Tube clamp
- > Latex, PVC, DEHP free



	<p>Description</p> <p>Product code: 7755643</p> <ul style="list-style-type: none"> > CH/FR 15 > ENFit connector > Wide lumen, gastric > External diameter: 4.8mm > Length: 35cm > Internal diameter: 3.6mm 	
<p>PERCUTANEOUS ENDOSCOPIC GASTROSTOMY USING THE PULL TECHNIQUE</p> <p>FREKA</p> <p>PEG STANDARD GASTRIC CH/FR 15</p> <p>For long-term intragastric feeding</p>	<p>Key Features</p> <ul style="list-style-type: none"> > Polyurethane tube with radio-opaque strip and silicone disc > External colour and size coded fixation plate (silicone) for securing tube (CH/FR 15 blue) > Puncture cannula with safety (air) valve > Introducer cone with double thread > Scalpel > ENFit connector > Interval markings every 2cm (from 2cm - 16cm) > Tube clamp > Latex, PVC, DEHP free 	



	<p>Description</p> <p>Product code: 7601365</p> <ul style="list-style-type: none"> > CH/FR 15 > 21cm shaft length > ENFit connector > Gastropexy device with suturing material > Trocar with peel away sheath > External diameter: 5.0mm 	
<p>DIRECT PUNCTURE GASTROSTOMY</p> <p>FREKA PEXACT II</p> <p>For medium to long-term intragastric feeding</p>	<p>Key Features</p> <ul style="list-style-type: none"> > Latex free silicone tube, balloon retained > Gastropexy device with suturing material, single handed operation > Trocar with peel away sheath, CH/FR 16 > Suturing material > Scalpel > Latex, PVC, DEHP free > Balloon with 90 day lifespan 	



JEJUNAL FEEDING TUBES



	<p>Description</p> <p>Product code: 7755644</p> <ul style="list-style-type: none"> > CH/FR 20 > ENFit connector > Wide lumen, gastric > External diameter: 6.6mm > Length: 35cm > Internal diameter: 5.0mm 	
<p>PERCUTANEOUS ENDOSCOPIC GASTROSTOMY USING THE PULL TECHNIQUE</p> <p>FREKA</p> <p>PEG STANDARD GASTRIC CH/FR 20</p> <p>For long-term intragastric feeding</p>	<p>Key Features</p> <ul style="list-style-type: none"> > Polyurethane tube with radio-opaque strip and silicone disc > External colour and size coded fixation plate (silicone) for securing tube (CH/FR 20 purple) > Puncture cannula with safety (air) valve > Introducer cone with double thread > Scalpel > ENFit connector > Interval markings every 2cm (from 2cm - 16cm) > Tube clamp > Latex, PVC, DEHP free 	

	<p>Description</p> <p>Product code: 7981832</p> <ul style="list-style-type: none"> > CH/FR 8, 270cm > ENFit connector > External diameter: 2.6mm > Internal diameter: 1.9mm 	
<p>NASOJEJUNAL FEEDING TUBE</p> <p>FREKA ENDOLUMINA</p>	<p>Key Features</p> <ul style="list-style-type: none"> > Tissue compatible tube made from polyurethane for endoscopic placement via the instrumental channel of the endoscope (cut to size required after placement) > Five staggered outlet ports with round closed end > ENFit connector with integral cap (attached after tube placement) > Radio-opaque > Interval markings every 1cm > Latex, PVC, DEHP free ^ <p>Re-routing catheter</p> <ul style="list-style-type: none"> > 35cm long, CH/FR 15 > Blue > Latex free <p>Accessory</p> <ul style="list-style-type: none"> > Freka Endo Introducer > 280cm long teflon coated stainless steel seldinger guide wire is available separately <p>Product code: 7980321</p>	

 <p>NASOJEJUNAL FEEDING TUBE</p> <p>FREKA</p> <p>EASY IN</p>	<p>Description Product code: 7981833</p> <ul style="list-style-type: none"> > Double lumen polyurethane tube suitable for intestinal feeding and gastric decompression/drainage > Intestinal: CH/FR 8, 270cm > ENFit connector > Internal diameter: 1.8mm > 5 lateral openings > Markings every 1cm > Gastric: CH/FR 16, 100cm > ENFit Connector > 4 lateral openings > Markings every 1cm 	
	<p>Key Features</p> <ul style="list-style-type: none"> > Tissue compatible polyurethane double lumen 270cm tube > Feeding lumen: 4 lateral openings, open ended, interval markings every 1cm > Gastric lumen: 5 large outlet ports, interval markings every 1cm > ENFit and funnel Y connector* > Pack contains 2 MCT oil ampoules for lubrication > Blue re-routing catheter > Latex, PVC, DEHP free ^ <p style="font-size: small;">* ISO80369 does not specify requirements for small bore connectors when used for gastric decompression (ISO/FDIS 80369-3:2016E)</p>	

 <p>FINE-NEEDLE CATHETER JEJUNOSTOMY FOR SURGICAL INSERTION</p> <p>FREKA</p> <p>SURGICAL JEJUNOSTOMY</p> <p>For postoperative enteral feeding following abdominal or laparoscopic procedures</p>	<p>Description Product code: 7755645</p> <ul style="list-style-type: none"> > CH/FR 9 > Length: 75cm > External diameter: 2.9mm > Internal diameter: 1.9mm > ENFit connector 	
	<p>Key Features</p> <ul style="list-style-type: none"> > Polyurethane tube with radio-opaque strip, open-ended ENFit connection with integral cap > External fixation plate (silicone) with clamp and suture holes (for securing to the abdomen with sutures) > Interval markings every 5cm > Two split cannulae (10.5cm and 13cm), blunt introducer in the longer split cannula > White clamp > Latex, PVC, DEHP free 	

 <p>NASOGASTRIC AND NASOJEJUNAL FEEDING TUBE</p> <p>FREKA</p> <p>TRELUMINA</p>	<p>Description Product code: 7981834</p> <ul style="list-style-type: none"> > Triple lumen polyurethane tube suitable for intestinal feeding and gastric decompression/drainage > Intestinal: CH/FR 9, ENFit, 150cm > Gastric: CH/FR 16, funnel connector, 95cm, pressure regulation lumen > Internal diameter: Feeding lumen 1.4mm, Aspiration lumen 4.1mm > Pressure regulation lumen 1.2mm 	
	<p>Key Features</p> <ul style="list-style-type: none"> > Tissue compatible polyurethane triple lumen 150cm tube > Feeding lumen: Flexible olive tip with 2 lateral outlets ports, open ended > Aspiration lumen: Five large outlet ports > Automatic pressure regulation lumen > Teflon coated guidewire 150cm > Radio-opaque strip along tube length > Interval markings every 1cm > Lumen identification marks > Pack contains 2 ampoules of MCT oil for lubrication of the tube and guidewire > Latex, PVC, DEHP free ^ 	

 <p>PERCUTANEOUS ENDOSCOPIC GASTROSTOMY USING THE PULL TECHNIQUE</p> <p>FREKA</p> <p>PEG UNIVERSAL INTESTINAL CH/FR 15</p> <p>For long-term intra-jejunal feeding</p>	<p>Description</p> <ul style="list-style-type: none"> > PEG Standard Gastric CH/FR 15, 35cm Product code: 7755643 > Intestinal CH/FR 9, 120cm Product code: 7755646 > Double lumen > External diameter: (CH/FR 9) 2.9mm > Internal diameter: (CH/FR 9) 1.9mm > ENFit connector with funnel connector for decompression 	
	<p>Key Features</p> <ul style="list-style-type: none"> > Tube and accessories (see Freka PEG universal gastric) > Additional radio-opaque polyurethane tube (120cm, CH/FR 15/ CH/FR 9) with teflon coated guidewire, olive tip and 5 lateral exit ports > The gastric tube acts as the introducer for the small intestine tube > Special connections available, depending on clinical requirements (straight or Y) > Interval markings <ul style="list-style-type: none"> • Gastric tube every 2cm • Intestinal tube every 5cm > Latex, PVC, DEHP free 	


ACCESSORY ITEMS

FREKA CLICK ADAPTOR FOR FREKA PEG INTESTINAL

DESCRIPTION

> For Freka PEG Intestinal CH/FR 9

Pack of 15
Product code: 7981389



FREKA Y CONNECTOR FOR FREKA PEG INTESTINAL

DESCRIPTION

> For Freka PEG Intestinal CH/FR 15

Pack of 5
Product code: 7981394




FREKA PEG CONNECTOR CH/FR 9 ENFIT

DESCRIPTION

> For Freka PEG CH/FR 9 (yellow)

Pack of 15
Product code: 7981385




FREKA PEG CONNECTOR CH/FR 15 ENFIT

DESCRIPTION

> For Freka PEG CH/FR 15 (blue)

Pack of 15
Product code: 7981386



FREKA STOMA MEASURING DEVICE

To be ordered separately when using the single set

DESCRIPTION

> Used to measure stoma tract length prior to button insertion

Pack of 1
Product code: 7750921




FREKA REPAIR SET FOR FREKA SURGICAL JEJUNOSTOMY, ENFIT

DESCRIPTION

> For Freka Surgical Jejunostomy CH/FR9

Pack of 1
Product code: 7981388




FREKA PEG CONNECTOR CH/FR 20 ENFIT

DESCRIPTION

> For Freka PEG CH/FR 20 (purple)

Pack of 15
Product code: 7981387




FREKA CAP ENFIT

DESCRIPTION

> ENFit syringe stopper

Pack of 30
Product code: 7981403



ACCESSORY ITEMS

FREKA FIXATION PLATE

DESCRIPTION

- > Freka PEG CH/FR 9

Pack of 5
Product code: 7903002




FREKA FIXATION PLATE

DESCRIPTION

- > Freka PEG CH/FR 15

Pack of 5
Product code: 7904002




FREKA FIXATION PLATE

DESCRIPTION

- > Freka PEG CH/FR 20

Pack of 5
Product code: 7751541




FREKA REPAIR SET

DESCRIPTION

- > For Freka PEG CH/FR 9/15/20
- > 1 external fixation plate
- > 1 clamp
- > 1 ENFit connector

Product	CH/FR 9	7981382
Code:	CH/FR 15	7981383
	CH/FR 20	7981384




FREKA EXTENSION SET (30CM)

DESCRIPTION

- > For use with balloon gastrostomies
- > Additional tube offering increased length
- > ENFit connector
- > Includes right-angle connector

Pack of 1
Product code: 7981397



FREKA EXTENSION SET (100CM)

DESCRIPTION

- > For use with balloon gastrostomies
- > Additional tube offering increased length
- > ENFit connector

Pack of 1
Product code: 7755693



CONE ADAPTOR

DESCRIPTION

- > Converts funnel end into ENFit end

Pack of 1 x15
Product code: 7751562




EASYBAG BOLUS ADAPTOR

DESCRIPTION

- > Resealable ENFit connection with ENPlus spiked end
- > For bolus feeding with EasyBag

Pack of 1 x 15
Product code: 7755694



ORDERING INFORMATION

GASTRIC TUBES AND BUTTONS

Item	Pack Size	Product Code
Freka Nasogastric Feeding Tube CH/FR 8, 80cm	1	7981839
Freka Nasogastric Feeding Tube CH/FR 8, 120cm	1	7981840
Freka Nasogastric Feeding Tube CH/FR 10, 120cm	1	7981842
Freka Nasogastric Feeding Tube CH/FR 12, 120cm	1	7981843
Freka GastroTube (CH/FR 15)	1	7755648
Freka Pexact II (CH/FR 15)	1	7601365
Freka Button Gastrostomy 1.1cm	1	7751241
Freka Button Gastrostomy 1.4cm	1	7751251
Freka Button Gastrostomy 1.7cm	1	7751261
Freka Button Gastrostomy 2.0cm	1	7751271
Freka Button Gastrostomy 2.5cm	1	7751281
Freka Button Gastrostomy 3.0cm	1	7751291
Freka Button Gastrostomy 3.5cm	1	7751301
Freka Button Gastrostomy 4.0cm	1	7751311
Freka Button Gastrostomy 4.5cm	1	7751321
Freka PEG standard gastric (CH/FR 9)	1	7755642
Freka PEG standard gastric (CH/FR 15)	1	7755643
Freka PEG standard gastric (CH/FR 20)	1	7755644

Freka Belly Button	CH/FR size/shaft length	CH/FR 12	CH/FR 14	CH/FR 16	CH/FR 18*	CH/FR 20*	CH/FR 24*
		Containing: 1 Freka Belly Button 1 Introducer (simple) 1 Syringe for balloon filling 10 ml Ampoule Aqua dest. 1 Lubricant 2 Gauze Pads	1.0 cm	7751612	7751622	7751632	7751642
	1.5 cm	7751613	7751623	7751633	7751643	7751653	7751663
	2.0 cm	7751614	7751624	7751634	7751644	7751654	7751664
	2.5 cm	7751615	7751625	7751635	7751645	7751655	7751665
	3.0 cm	7751616	7751626	7751636	7751646	7751656	7751666
	3.5 cm	7751617	7751627	7751637	7751647	7751657	7751667
	4.0 cm	7751618	7751628	7751638	7751648	7751658	7751668
	4.5 cm	7751619	7751629	7751639	7751649	7751659	7751669
	5.0 cm	7751620	7751630	7751640	7751650	7751660	7751670

* Freka Belly Button sizes CH/FR 18 CH/FR 20 CH/FR 24 may have longer delivery times.

JEJUNAL FEEDING TUBES

Item	Pack Size	Product Code
Freka Endolumina	1	7981832
Freka Easy In	1	7981833
Freka Trelumina	1	7981834
Freka Surgical Jejunostomy	1	7755645
Freka Intestinal Tube (CH/FR 9) to be used with Freka PEG CH/FR 15	1	7755646

ACCESSORY ITEMS

Item	Pack Size	Product Code
Freka Endo Introducer (280cm)	1	7980321
Freka Click adaptor for Freka PEG intestinal CH/FR 9	1 x 15	7981389
Freka Y Connector for Freka Intestinal CH/FR 15	1 x 5	7981394
Freka Stoma Measuring Device	1	7750921
Freka PEG connector CH/FR 9 ENFit	1 x 15	7981385
Freka PEG connector CH/FR 15 ENFit	1 x 15	7981386
Freka PEG connector CH/FR 20 ENFit	1 x 15	7981387
Fixation plate for PEG CH/FR 9	1 x 5	7903002
Fixation plate for PEG CH/FR 15	1 x 5	7904002
Fixation plate for PEG CH/FR 20	1 x 5	7751541
Freka Cap ENFit	1 x 30	7981403
Freka Repair Set for Freka Surgical Jejunostomy, ENFit	1 x 15	7981388
Repair set for Freka PEG CH/FR 9	1	7981382
Repair set for Freka PEG CH/FR 15	1	7981383
Repair set for Freka PEG CH/FR 20	1	7981384
Freka Extension set (30cm) (Right Angled)	1 x 1	7981397
Freka Extension set (100cm)	1 x 1	7755693
Cone adaptor	1 x 15	7751562
EasyBag Bolus adaptor	1 x 15	7755694
ENFit syringe stopper	1 x 30	7981403

^ Please note that whilst our Freka Endolumina, Trelumina and Easy In tubes are all DEHP free, the packs contain re-routing catheters/tubes which are disposed of upon tube placement, made from PVC with DEHP.



**FRESENIUS
KABI**

caring for life

Fresenius Kabi Limited,
Cestrian Court, Eastgate Way,
Manor Park, Runcorn,
Cheshire WA7 1NT
tel: 01928 533533 fax: 01928 533520
email: scientific.affairsUK@fresenius-kabi.com
www.fresenius-kabi.co.uk

Date of preparation: October 2019 EN00233

- ® The signs/names marked with ® are registered trademarks of the Fresenius Group in selected countries.
- ® Freka is a registered trademark of Fresenius Kabi AG. Fresenius Kabi Ltd is an authorised user.
- ® ENFit is a registered trademark of GEDSA, Inc.
- © Fresenius Kabi Ltd. June 2017